Analysis Platform based on Django & Keras

7106056041 靳開翔




Outline

- **研究動機**
- ■研究目標
- ■資料初探
- ■前處理
- 實作過程
- ■實驗結果



研究動機





研究目標

- 1. 使用兩種以上的演算法,對ECG訊號進行心臟疾病辨識。
 - 2. 能夠辨識出

 - 心房期外收縮心跳 (Atrial premature contraction · APC)

 - 左分枝束阻斷心跳(Left bundle branch block beat, LBBB)
 - 心室跳脫心跳 (Ventricular escape beat)
 - 心室撲動 (Ventricular flutter wave) 等疾病。
 - 3. 使用Django架置這個有EA系統的網站



MIT-BIH ARRHYTHMIA DATABASE

Follow the links in this table to go directly to any record:

<u>100 104 108 113 117 122 201 207 212 217 222 231</u>

101 105 109 114 118 123 202 208 213 219 223 232

102 106 111 115 119 124 203 209 214 220 228 233

103 107 112 116 121 200 205 210 215 221 230 234

Record 100 (MLII, V5; male, age 69)

Medications: Aldomet, Inderal

Beats	Before 5:00	After 5:00	Total
Normal	367	1872	2239
APC	4	29	33
PVC	-	1	1
Total	371	1902	2273



資料初探

■ 100.atr


100.dat

100.hea

Time	Sample #	Туре	Sub	Chan	Num	Symbol	
0:00.050	18	+	0	0	0	· or N	Normal beat
0:00.214	77	N	0	0	0	L	Left bundle branch block beat
0:01.028	370	N	0	0	0	R	Right bundle branch block beat
0:01.839	662	N	0	0	0	A	Atrial premature beat
0:02.628	946	N	0	0	0	a	Aberrated atrial premature beat
0:03.419	1231	N	0	0	0	J	Nodal (junctional) premature beat
0:04.208	1515	N	0	0	0	S	Supraventricular premature beat
0:05.025	1809	N	0	0	0	V	Premature ventricular contraction
						F	Fusion of ventricular and normal beat
0:05.678	2044	A	0	0	0	[Start of ventricular flutter/fibrillation
0:06.672	2402	N	0	0	0	!	Ventricular flutter wave
0:07.517	2706	N	0	0	0]	End of ventricular flutter/fibrillation
0:08.328	2998	N	0	0	0	е	Atrial escape beat
0:09.117	3282	N	0	0	0	j	Nodal (junctional) escape beat
0:09.889	3560	N	0	0	0	E	Ventricular escape beat
0:10.728	3862	N	0	0	0	/	Paced beat
0:11.583	4170	N	0	0	0	f	Fusion of paced and normal beat
0:12.406	4466	N	0	0	0	х	Non-conducted P-wave (blocked APB)
0:13.233	4764	N	0	0	0	Q	Unclassifiable beat
0:14.056	5060	N	0	0	0		Isolated QRS-like artifact



前處理



So and Chan

	id	Pamp	Ramp	Tamp	PR	QRS	QT	RRI	dRRI	label
0	1	-0.240	0.810	-0.270	133.333	47.222	352.778	788.889	22.222	N
1	2	-0.230	0.820	-0.305	183.333	41.667	344.444	791.667	2.778	N
2	3	-0.195	0.885	-0.280	138.889	44.444	372.222	788.889	2.778	N
3	4	-0.200	0.945	-0.270	138.889	44.444	366.667	816.667	27.778	N
4	5	-0.205	0.845	-0.305	163.889	41.667	369.444	652.778	163.889	Α



MinMaxScaler:


Feature scaling is used to bring all values into the range [0,1].


$$X' = rac{X - X_{\min}}{X_{\max} - X_{\min}}$$



LSTM Parameters

Reference: ECG-BASED BIOMETRICS USING RECURRENT NEURAL NETWORKS







分類結果

Classifier	Naïve Bayes	Linear SVC	KNN	Random Forest	LSTM
Accuracy	48%	45%	89%	96%	97%



Django 架構


新架構





Django 操作介面

ECG Analysis Platform



© Jin 2018



Django 操作介面

